Italian  
Medical  
Touch

**FLO** ————— **MIKO**

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**MIKO** is an electrical fed compressor with aerosol therapy atomiser.

The instrument is designed for easy transport and holding and is recommended for atomising antibiotics and bronchodilator drugs.

The high thermal insulated plastic body complies with European Safety Rules.

The medical device is designed for continuous use.

## GENERAL WARNING

**READ INSTRUCTION MANUAL CAREFULLY BEFORE USE**

**DRUG ADMINISTRATION MUST BE UNDER MEDICAL CONTROL**

**THE INSTRUMENT MUST NOT BE DISASSEMBLED  
FOR A TECHNICAL SERVICE ALWAYS CONTACT CA-MI**

## IMPORTANT SAFETY RULES

1. On opening the packaging, check the integrity of the appliance, paying particular attention to the presence of damage to the plastic parts, which may make access possible to internal live parts and also to breakage and / or peeling of the power supply cable. **In these cases don't connect the plug to the electric socket. Carry out these controls before each use;**
2. before connecting the appliance always check that the electric data indicated on the data label and the type of plug used, correspond to those of the mains electricity to which it's to be connected;
3. If the plug supplied with the appliance is incompatible with the mains electricity socket, contact qualified staff for replacement of the plug with a suitable type. The use of simple or multiple and / or extension adapters is not generally recommended.  
Whenever their use is indispensable, use those in compliance with safety regulations, however paying attention not to exceed the maximum power supply limits, which are indicated on the adapters and extensions;
4. Never leave the appliance inserted if not necessary disconnect the plug from the mains power supply when it is not being used;
5. Respect the safety regulations indicated for electrical appliances and particularly:
  - Only use original accessories and components;
  - Never immerge the appliance into water;
  - Position the appliance on flat stable surfaces;
  - Position the device in a way that the air inlets on the back aren't obstructed;
  - Never use the device in environments which have anaesthetic mixtures inflammable with air, oxygen or nitric oxide;
  - Don't touch the device with wet hands and always prevent the appliance coming into contact with liquids;
  - The use of this device by children and / or incompetent person always requires the careful surveillance of an adult in possession of their full mental faculties;
  - The medical device, and most of all the nebulae, must be kept out of children's reach as it contains small parts that could be swallowed;
  - Don't leave the appliance connected to the power supply socket when not in use;
  - Don't pull the power supply cable to disconnect the plug remove the plug from the mains socket correctly;
  - Preserve and use the medical device in environments protected from atmospheric factors and at a distance from heat sources;
6. For repairs, exclusively contact CA-MI technical service and request the use of original spare parts. Failure to comply with the above can jeopardise the safety of the device;
7. **This medical device must be destined exclusively for the use for which it has been designed and described in this manual. It must therefore be used as an aerosol therapy system.**  
Any different use must be considered incorrect and therefore dangerous; the manufacturer cannot be considered liable for damage caused by improper, incorrect and / or unreasonable use or if the appliance is used in electrical plants that are not in compliance with the regulations in force;

8. Particular precautions must be made concerning electromagnetic compatibility. The medical device must be installed and used according to information supplied with the accompanying documents;
9. Some components of the device are small enough to be swallowed by children: therefore keep the device out of children's reach;
10. Remember to:
  - Only use this device with medicines prescribed by your doctor;
  - Carry out the treatment only using the accessory indicated by the doctor according to the pathology.

**REMOVE THE DEVICE FROM CARRYING BAG BEFORE USE**

**INSTRUCTIONS FOR THE CORRECT DISPOSAL OF THE PRODUCT FOR THE PURPOSE AND EFFECT OF EUROPEAN COMUNITY DIRECTIVE 2002/96/EC::**

At the end of its useful life the product must not be disposed of with domestic waste. It may be delivered to the relevant separate waste collection centres prepared by municipal administrations, or at retailers that supply this service. Disposing of the product separately prevents possible negative consequences for the environment and for health, deriving from inadequate disposal. It also allows the recovery of materials of which it's made up in order to obtain an important saving of energy and resources. In order to point out the obligation to dispose of the electromedical appliances separately, the product has the sign of the crossed mobile waste container applied to it.


**RULES FOR RETURNING AND REPAIRING**

**COMPLYING WITH THE NEW EUROPEAN RULES, CA-MI INDICATES THE IMPORTANT POINTS TO PROTECT INSTRUMENT AND OPERATORS HYGIENE. THESE RULES MUST BE RESPECTED IN ORDER TO GUARANTEE HYGIENE AND SAFETY TO ALL THE PEOPLE OPERATING WITH THE INSTRUMENT TO OBTAIN QUALITY AND WELL BEING.**

Every returned instrument will be hygienically checked before repairing. If CA-MI finds instrument not suitable for repairing due to clear signs of internal or external contamination, the same will be returned to customer with specification of NOT REPAIRED INSTRUMENT, accompanied by an explanation letter. CA-MI will decide if contamination is due to bad functioning or misuse. If contamination is due to bad functioning, CA-MI will substitute the instrument, only if SALE RECEIPT and STAMPED GUARANTEE accompany the same. CA-MI is not responsible for contaminated accessories, they will be substitute at customer's expenses. For this reason it is **COMPULSORY** to carefully disinfect the external part of the instrument and accessories with a cloth soaked in methylated spirits or hypochlorite-based solutions. Put the instrument and accessories in a bag with indication of disinfecting. We also request to specify the kind of fault, in order to speed up repairing procedures. To this end, please read the instructions carefully in order to avoid damaging the equipment through improper use. Always specify the fault encountered so that CA-MI can establish whether it falls into the category of the faults covered by the guarantee.


## TECHNICAL CHARACTERISTICS

TYPOLOGY (MDD 93/42/EEC)	Class IIa Medical Device
MODEL	<b>MIKO</b>
POWER FEEDING	230V~ / 50Hz
POWER CONSUMPTION	170VA
FUSE	F 1 x 1.6A 250V
MAX PRESSURE	220 kPa (2.2 Bar)
MAX AIR FLOW	13 l/min
OPERATING PRESSURE	95 kPa (0.98 Bar)
OPERATING AIR FLOW	4.8 l/min a 95kPa
NEB-RATE (with 4ml of 0.9% NaCl solution)	0.35 ml/min
MMAD	2.44 $\mu$ m
GSD	2.87
WEIGHT	1.45 Kg
SIZE	148 x 171 x 114 (h)
NOISE LEVEL (measured as specifications of EN 3544-1)	Approx. 56dB (A)
DUTY CYCLE (to 40°C and 110% operating voltage)	NOT-STOP OPERATED
MIN. CAPACITY NEBULIZER	2ml
MAX. CAPACITY NEBULIZER	6ml
WORKING CONDITION	Room temperature: 10 $\div$ 40°C Room humidity percentage: 20 $\div$ 85% RH Altitude: 0 $\div$ 2000m s.l.m.
CONSERVATION CONDITION AND TRASPORT	Room temperature: -25 $\div$ 70°C Room humidity percentage: 10 $\div$ 95% RH

MMAD = Mass Median Aerodynamic Diameter  
GSD = Geometric Standard Deviation

## SYMBOLS

	Class II isolation equipment
<b>CE 0123</b>	CE marking in conformity with EC directive 93/42/EEC
	Warning, consult the instruction manual
	To Preserve in place coolness and dry land
	Conservation temperature: -25 ÷ 70°C
	Type B equipment
	Fuse
~	Alternate Current
Hz	Mains Frequency
I	ON
O	OFF

## CLEANING DEVICE

Use a soft dry cloth with not – abrasive and not – solvent detergents.

**MAKE SURE THAT INTERNAL PARTS OF THE APPLIANCE DON'T COME INTO CONTACT WITH LIQUIDS AND THE PLUG IS NOT INSERTED**

## CELANING THE ACCESSORIES

**Procced as follows to disinfect the accessories:**

1. Turn the upper part of the nebulizer in an anti-clockwise direction;
2. Disconnect the internal pisper at the base of the nebulizer using the fingers;
3. Use denatured alcohol or a hypochlorite - based solution, easily found at chemist, to clean the accessories;
4. After using the appliance, disassemble the nebulizer and clean all parts in warm water, rinse carefully and remove excess water using a soft cloth and leave to dry in a clean place.

**DO NOT BOIL OR PUT IN AUTOCLAVE THE ACCESSORIES**

After each treatment clean thoroughly each component of the nebulizer (except air tube) removing medication residual and possible impurities. Clean all parts in warm water.  
Rinse thoroughly making sure that all deposits are washed away and let dry.

<b>Guidance and manufacturer's declaration – Electromagnetic Emissions</b> (as request by regulation EN 60601-1-2:2001)		
The <b>MIKO</b> aerosol is intended for use in the electromagnetic environment specified below. The customers or the user of the <b>MIKO</b> aerosol should assure that it's used in such an environment.		
<b>Emissions Test</b>	<b>Compliance</b>	<b>Electromagnetic environment - guidance</b>
Irradiated / Conducted emissions CISPR11	Group 1	The MIKO aerosol only used RF energy only for its internal functioning. Therefore its RF emissions are very low and are not cause interference in proximity of any Electronic appliances.
Irradiated / Conducted emissions CISPR11	Class [B]	The MIKO aerosol can be used in all environments, including domestic and those connected directly to the public mains distribution that supplies power to environments used for domestic scopes.
Harmonic emissions IEC/EN 61000-3-2	Class [A]	
Voltage fluctuations / flicker emissions IEC/EN 61000-3-3	Complies	

<b>Guidance and manufacturer's declaration – Electromagnetic Emissions</b> (as request by regulation EN 60601-1-2:2001)		
The <b>MIKO</b> aerosol is intended for use in the electromagnetic environment specified below. The customers or the user of the <b>MIKO</b> aerosol should assure that it's used in such an environment.		
<b>Immunity Test</b>	<b>Compliance</b>	<b>Electromagnetic environments - guidance</b>
Electrostatic discharge (ESD) IEC/EN 61000-4-2	$\pm 6\text{kV}$ on contact $\pm 8\text{kV}$ in air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient / burst IEC/EN 61000-4-4	$\pm 2\text{kV}$ power supply	Mains power quality should be that of a typical commercial environment or hospital
Surge IEC/EN 61000-4-5	$\pm 1\text{kV}$ differential mode	Mains power quality should be that of a typical commercial environment or hospital
Loss of voltage, brief voltage interruptions and variations IEC/EN 61000-4-11	$5\%U_T$ for 0.5 cycle $40\%U_T$ for 05 cycle $70\%U_T$ for 25 cycle $<5\%U_T$ for 5 sec	Mains power quality should be that of a typical commercial environment or hospital If the user of the MIKO aerosol request that the appliance operates continuously, the use of a continuity unit is recommended.
Magnetic field IEC/EN 61000-4-8	3A/m	The power frequency magnetic field should be measured in the intended installation location to assure that it's sufficiently low.
Conducted Immunity IEC/EN 61000-4-6	3Vrms 150kHz to 80MHz (for appliances that aren't life - supporting)	-
Irradiated Conducted IEC/EN 61000-4-3	3V/m 80MHz to 2.5 GHz (for appliances that aren't life - equipment)	-
Note $U_T$ is the value of the power supply voltage		

## INSTRUCTION FOR USE


- Extract the power supply cable and insert the plug into the mains socket. It is recommended to unwind the entire length of the power supply cable to prevent dangerous overheating. If the power supply cable is damaged and must be replaced contact the CA-MI technical service;
- Open the nebulizer **2** by unscrewing the lid;
- Pour the medicine prescribed by the doctor into the nebulizer;
- Re-close the nebulizer, re-screwing the lid;
- Connect air pipe **5** to the air exit well **4**;
- Connect the other end of the pipe to the connection in the lower part of the nebulizer;
- Connect the desired accessory to the nebulizer: child mask or adult mask or mouth-piece;
- Ensure that the supplied air filter **(6)** is present;
- Press switch **1** on position **I** to proceed with nebulization;
- On completing of nebulization, press the switch on position **0** and remove the plug from the socket;
- Wash the nebulizer and its accessories as indicated in the cleaning charter;.
- Place the cable and accessories inside the box.

Always use the nebulizer facing upwards so that substances and / or medicines cannot escape from the nebulizer during the normal use.

**NEVER INHALE IN HORIZONTAL POSITION  
NEVER BEND THE NEBULIZER OVER 60°**

- 1 – Air Tube
- 2 – Nebulizer Tank
- 3 – Nebulization Nozzle
- 4 – Nebulizer Top
- 5 – Mouthpiece
- 6 – Adult Mask
- 7 – Pediatric Mask

## STANDARD ACCESSORIES

ACCESSORIES
HI-FLO KIT (Nebulizer HI-FLO, Adult Mask, Pediatric Mask, Air Tube and Mouth-piece)
Air Filter (+n°1 spare part)
Carrying Bag (on request)

For each individual patient it's recommended to use the nebulizer for 6 months or for a maximum of 120 treatments. The nebulizer must be replaced after a long period of inactivity, if it is deformed or broken, or if the nebulizer nozzle is blocked by dry medicine, dust, ecc..

**Only use the original nebulizer supplied by CA-MI with the device**

The air filter must be replaced every 25 hours of functioning or when it result particulary worn.

For replacement, lift the filter and replace with a new one.

Only use original CA-MI filter

**DON'T USE THE DEVICE WITHOUT AIR FILTER**

## MAINTENANCE

The **MIKO** atomiser does not need maintenance or lubrication.

It is necessary to check functioning and instrument before every use. Unpack the instrument and **always check** integrity of plastic parts and feeding cable., they might have been damaged during previous use.

Connect cable to electrical network and turn switch on.

Close the compressor mouthpiece with a finger and verify that loud noises are not present, these can indicate wrong functioning. Verify that the atomiser is not damaged by previous use (it was badly put away or badly knocked). A protection fuse (**F 1.6 A 250V**) not reachable from exterior protects the instrument.

For fuse replacing, please make reference to manufacturer technical personnel.

Fault type	Cause	Solution
1. Low Nebulization	Clogged Nebulizer Tank	Clean and disinfect the nebulizer tank as explained in the instruction manual
2. Low Nebulization	Clogged Nebulizer Tank	If cleaning was not succesful change cruet
3. Absence of Nebulization	Clogged Nebulizer Tank	Check that the nebulizer contains medication; Make sure that the nebulizer is not clogged; Check the connection between the compressor air outlet port and the accessories
4. Slow Nebulization	Highly drug	Dilute drug in physiological liquid
5. Noisy Device	Extended use	Call retainer or manufacturer CA-MI
<b>Fault 1 - 2 - 3 - 4 - 5</b>	<b>No solution with previous items</b>	<b>Call retainer or manufacturer CA-MI</b>

If the unit does not nebulizer once the above conditions have been checked, we suggest to contact your dealer or technical service CA-MI.

**BEFORE EVERY CHECKING OPERATION, IN CASE OF ANOMALIES OR BAD FUNCTIONING,  
PLEASE CONTACT CA-MI TECHNICAL SERVICE  
CA-MI DOES NOT GIVE GUARANTEE IF INSTRUMENT, AFTER THE TECHNICAL SERVICE  
CHECKING, APPEARS TO BE TAMPERED**

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### WARRANTY

CA-MI warrants its products for **36 months** after purchasing date. In front of this warranty, CA-MI will be obliged only to repair or substitute free of charge the products or parts of them that, after verification effected on our factory, or our authorized Service Center, by the Technical Service, results defective. The product must be accompanied by a description of the defect. The warranty, with exclusion of responsibility for direct and indirect damages, it is thought limited to the solos defects of material or workmanship and it stops having effect when the device results however gotten off, tampered or sheltered out of the Factory or from the Authorized Service center.

The commodity always travels to risk and danger of the buyer, without any responsibility of CA-MI for damages caused by the transport or dismay from the vector.

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
**TIMRO E FIRMA RIVENDITORE  
AUTHORIZED DEALER STAMP AND  
SIGNATURE**

**FOR WARRANTY REPARATIONS THE  
POSTCARD WILL BE RETURNED WITH  
THE DEVICE**

**CA.MI. di Attolini Mario & C. s.n.c.**  
Via Ugo La Malfa nr.31  
43010 Pilastro (PR) Italia  
Tel. +39 0521 / 637133 – 631138  
Fax. +39 0521 / 639041  
E-mail: [vendite@ca-mi.it](mailto:vendite@ca-mi.it)  
[export@ca-mi.it](mailto:export@ca-mi.it)